



**10X MOPS Buffer for EZ Rich
Defined *E. coli* Growth Medium
(Cat. No. D-0508-10)**

FOR RESEARCH USE ONLY

BACKGROUND

MOPS- 3-(N-morpholino) propanesulfonic acid. MOPS ($pK_a = 7.20$) is an excellent buffer for many biological systems at near-neutral pH, MOPS Minimal is designed for culturing *E. coli* and works very well with our MDST™ 42 strains.

REAGENTS	VOLUME
----------	--------

10X MOPS	100ml
0.132M K_2PO_4	10ml

MEDIA PREPARATION

Add reagents plus carbon source to approximately 800 ml pure water; stir to completely mix. Adjust volume to 1L. If measured pH is not 7.2, adjust by dropwise addition of 16M HCl. Filter sterilize using 0.2 μ m vacuum filter into sterile container.

PLATE PREPARATION

Add 13g/L Agar to approximately 880ml water (this volume will vary depending on the volume of carbon source(s) added), autoclave for 35 min at 120°C. While agar solution is autoclaving, combine the MOPS reagents. If measured pH is not 7.2, adjust by dropwise addition of 16M HCl. Filter sterilize using a 0.2 μ m vacuum filter. Once the agar solution has cooled to 50°C add the MOPS solution, mix aseptically, and pour plates.

STORAGE CONDITIONS

Store 10X MOPS reagent at -20 °C.
Store 0.132 K_2PO_4 at room temperature.